# Bioreset Smart



The system can be placed directly inside the environment under treatment. To easily move it in different areas, it can be equipped with a trolley that includes a balance. It works on stand-alone mode or can be managed through a remote connection via ModBus TCP/IP.

### **FEATURES**

- Ergonomic Italian design for maximum operator
- Integrated 7" colour touch panel, simple and intuitive
- Operating stand-alone
- Communication via Ethernet connection (ModBus TCP/IP) for remote control
- · Software with two account levels, recipes and reports storage
- Adjustable blower up to 200 m<sup>3</sup>/h
- · Reports exported in pdf format via USB port

# **TYPICAL APPLICATIONS**

- · Clean room for sterile production
- BSL-3 Laboratory



# Bioreset Smart

#### **OPERATIONS**

- Set-up the Bioreset Smart inside the area or the application to be decontaminated
- Turn off air conditioning / HVAC system
- Secure area with warning signs for personal safety
- Start the recipe on Bioreset Smart and leave it inside the area until the cycle automatically ends
- Aerate the room/application using the air conditioning / HVAC system
- For personnel safety, upon cycle end check H<sub>2</sub>O<sub>2</sub> gas residue before anybody can enter the room
- Remote control allows to start and stop the Bioreset Smart and to receive alarm signals for different failures.


# **OPTIONAL ACCESSORIES**

- · Catalyst for exhaust or return air
- Special pipe dimensions on request
- Safety devices for personnel
- Cables for remote connection
- T/RH% probe package for automatic cycle
- Movable trolley with built-in balance and safety support for H<sub>2</sub>O<sub>2</sub> bottles of 1 and 5 litres
- Turboflow, a versatile system for ultra-fast H<sub>2</sub>O<sub>2</sub> distribution



| Technical data                      |  |
|-------------------------------------|--|
| Max treatable volume*               | Up to 500 m <sup>3</sup>                                       |
| Nominal blower airflow              | Up to 200 m³/h, adjustable                                     |
| T/RH% probe                         | Optional<br>0% 100% RH – step 0,1%<br>-40°C +60°C – step 0,1°C |
| Balance                             | Optional on trolley<br>0 g 10.000 g – step 1 g                 |
| Peristaltic pump                    | Integrated<br>3 12 g/min – step 0,5                            |
| Remote control                      | PC or PLC  |
| Communication protocol              | ModBus TCP/IP  |
| Ø in / out connection               | 60 mm or tri-clamp 2"1/2                                       |
| Materials                           | Stainless steel AISI 316 and ABS                               |
| Power supply                        | 230V / 50Hz or 110V / 60Hz                                     |
| Max consumption                     | 1 kW   |
| Dimensions (L x P x H)              | 325 x 570 x 350 mm   |
| Weight                              | 24 kg  |
| Dimensions with trolley (L x P x H) | 485 x 700 x 1.200 mm   |
| Weight with trolley                 | 74 kg  |

<sup>\*</sup> clean, dry, sealed enclosures



### Legenda

- 1 USB interface
- 2 Ethernet interface
- **3** H<sub>2</sub>O<sub>2</sub> inlet
- 4 Power supply
- 5 Balance interface
- 6 T/RH% probe
- 7 Air oulet
- 8 Touch panel
- 9 Air inlet

